OPTIKA

M I C R O S C O P E S

I T A L Y

Ver. 3.0.0

B-500DK ***B-500DK-R***

OPERATION MANUAL

MANUAL DE INSTRUCCIONES

OPTIKA MICROSCOPES - ITALY

www.optikamicroscopes.com - info@optikamicroscopes.com

1.0 B-500DK DARKFIELD MICROSCOPE	page 3
2.0 PRINCIPLES OF IMMERSION MICROSCOPY	page 3
3.0 PRINCIPLES OF DARKFIELD ILLUMINATION	page 6
3.1 Darkfield Microscopy at High Magnifications	
4.0 B-500DK DARKFIELD MICROSCOPE CONFIGURATION	page 10
4.1 Troubleshooting	
5.0 TECHNICAL TIPS FOR OIL IMMERSION MICROSCOPY	page 20
6.0 ELECTRICS	page 21
7.0 USING THE RECHARGEABLE BATTERIES	page 22
8.0 RECYCLING AND RECOVERY	page 26

B-500DK is a darkfield system specific for blood analysis with a 1.36 - 1.25 N.A. special extra efficient darkfield condenser and a 100X plan-achromatic objective with adjustable iris diaphragm. The X-LED illumination ensures the high level of light intensity typically needed in high magnification darkfield techniques.

In order to correctly use this microscope, one has to gain some familiarity with a) oil immersion technique and b) darkfield technique.

In the following manual we present the basics of these methods (chapters 2 and 3) and then we give a step-by-step guide to configuration of B-500DK (chapter 4). General tips for immersion microscopy are also given (chapter 5).

2.0 PRINCIPLES OF IMMERSION MICROSCOPY

The ability of a microscope objective to capture deviated light rays from a specimen is dependent upon both the numerical aperture and the medium through which the light travels.

An objective's numerical aperture is directly proportional to the refractive index of the imaging medium between the coverslip and the front lens, and also to the sin of one-half the angular aperture of the objective. Because sin cannot be greater than 90 degrees, the maximum possible numerical aperture is determined by the refractive index of the immersion medium. Most microscope objectives use air as the medium through which light rays must pass between the coverslip protecting the sample and front lens of the objective. Objectives of this type are referred to as dry objectives because they are used without liquid imaging media. Air has a refractive index of 1.0003, very close to that of a vacuum and considerably lower than most liquids, including water ($n = 1.33$), glycerin ($n = 1.470$) and common microscope immersion oils (average $n = 1.515$). In practice, the maximum numerical aperture of a dry objective system is limited to 0.95, and greater values can only be achieved using optics designed for immersion media.


Figure 1 Oil immersion and numerical aperture

The principle of oil immersion is demonstrated in Figure 1 where individual light rays are traced through the specimen and either pass into the objective or are refracted in other directions. Figure 1 (a) illustrates the case of a dry objective with five rays (labeled 1 through 5) shown passing through a sample that is covered with a coverslip. These rays are refracted at the coverslip-air interface and only the two rays closest to the optical axis (rays 1 and 2) of the microscope have the appropriate angle to enter the objective front lens. The third ray is refracted at an angle of about 30 degrees to the coverslip and does not enter the objective. The last two rays (4 and 5) are internally reflected back through the coverslip and, along with the third ray, contribute to internal reflections of light at glass surfaces that tend to degrade image resolution. When air is replaced by oil of the same refractive index as glass, shown in Figure 1(b), the light rays now pass straight through the glass-oil interface without deviation due to refraction. The numerical aperture is thus increased by the factor of n , the refractive index of oil.

Microscope objectives designed for use with immersion oil have a number of advantages over those that are used dry. Immersion objectives are typically of higher correction (either fluorite or apochromatic) and can have working numerical apertures up to 1.40 when used with immersion oil having the proper dispersion and viscosity. These objectives allow the substage condenser diaphragm to be opened to a greater degree, thus extending the illumination of the specimen and taking advantage of the increased numerical aperture.

A factor that is commonly overlooked when using oil immersion objectives of increased numerical aperture is limitations placed on the system by the substage condenser. In a situation where an oil objective of $NA = 1.40$ is being used to image a specimen with a substage condenser of smaller numerical aperture (1.0 for example), the lower numerical aperture of the condenser overrides that of the objective and the total NA of the system is limited to 1.0, the numerical aperture of the condenser.

Modern substage condensers often have a high degree of correction with numerical aperture values ranging between 1.0 and 1.40. In order to effectively utilize all the benefits of oil immersion, the interface between the substage condenser front lens and the underside of the microscope slide containing the specimen should be also be immersed in oil. An ideal system is schematically diagramed in Figure 2, where immersion oil has been placed at the interfaces between the objective front lens and the specimen slide and also between the front lens of the condenser and the underside of the specimen slide.

This system has been termed a Homogeneous Immersion System and it is the ideal situation to achieve maximum numerical aperture and resolution in an optical microscope.


Figure 2 Homogeneous immersion system

In this case, the refractive index and dispersion of the objective front lens, immersion oil, substage condenser front lens, and the mounting medium are equal or very near equal.

In this ideal system, an oblique light ray can pass through the condenser lens and completely through the microscope slide, immersion oil, and mounting medium undeviated by refraction at oil-glass or mounting medium-glass interfaces.

When using high-power achromat oil immersion objectives, it is sometimes permissible to omit the step of oiling the condenser top lens. This is because the condenser aperture diaphragm must often be reduced with lesser-corrected objectives to eliminate artifacts and provide optimum imaging. The reduction in diaphragm size reduces the potential increase in numerical aperture (provided by oiling the condenser lens) so the loss in image quality under these conditions is usually negligible.

Darkfield microscopy is a specialized illumination technique that capitalizes on oblique illumination to enhance contrast in specimens that are not imaged well under normal brightfield illumination conditions.

All of us are quite familiar with the appearance and visibility of stars on a dark night, this despite their enormous distances from the earth. Stars can be seen because of the stark contrast between their faint light and the black sky.

This principle is applied in darkfield (also called darkground) microscopy, a simple and popular method for making unstained objects clearly visible. Such objects often have refractive indices very close in value to that of their surroundings and are difficult to image in conventional brightfield microscopy. For instance, many small aquatic organisms have a refractive index ranging from 1.2 to 1.4, resulting in a negligible optical difference from the surrounding aqueous medium. These are ideal candidates for darkfield illumination.

Darkfield illumination requires blocking out of the central light which ordinarily passes through and around (surrounding) the specimen, allowing only oblique rays from every azimuth to “strike” the specimen mounted on the microscope slide. The top lens of a simple Abbe darkfield condenser is spherically concave, allowing light rays emerging from the surface in all azimuths to form an inverted hollow cone of light with an apex centered in the specimen plane. If no specimen is present and the numerical aperture of the condenser is greater than that of the objective, the oblique rays cross and all such rays will miss entering the objective because of their obliquity. The field of view will appear dark.

The darkfield condenser/objective pair illustrated in Figure 3 is a high-numerical aperture arrangement that represents darkfield microscopy in its most sophisticated configuration, which will be discussed in detail below. The objective contains an internal iris diaphragm that serves to reduce the numerical aperture of the objective to a value below that of the inverted hollow light cone emitted by the condenser. The cardioid condenser is a reflecting darkfield design that relies on internal mirrors to project an aberration-free cone of light onto the specimen plane.

When a specimen is placed on the slide, especially an unstained, non-light absorbing specimen, the oblique rays cross the specimen and are diffracted, reflected, and/or refracted by optical discontinuities (such as the cell membrane, nucleus, and internal organelles) allowing these faint rays to enter the objective. The specimen can then be seen bright on an otherwise black background. In terms of Fourier optics, darkfield illumination removes the zeroth order (unscattered light) from the diffraction pattern formed at the rear focal plane of the objective. This results in an image formed exclusively from higher order diffraction intensities scattered by the specimen.


Figure 3 Cardioid darkfield condenser

Ideal candidates for darkfield illumination include minute living aquatic organisms, diatoms, small insects, bone, fibers, hair, unstained bacteria, yeast, and protozoa. Non-biological specimens include mineral and chemical crystals, colloidal particles, dust-count specimens, and thin sections of polymers and ceramics containing small inclusions, porosity differences, or refractive index gradients. Care should be taken when preparing specimens for darkfield microscopy because features that lie above and below the plane of focus can also scatter light and contribute to image degradation. Specimen thickness and microscope slide thickness are also very important and, in general, a thin specimen is desirable to eliminate the possibility of diffraction artifacts that can interfere with image formation.

For more precise work and blacker backgrounds, you may choose a condenser designed especially for darkfield, i.e. to transmit only oblique rays. There are several varieties: “dry” darkfield condensers with air between the top of the condenser and the underside of the slide—and immersion darkfield condensers which require the use of a drop of immersion oil (some are designed to use water instead) establishing contact between the top of the condenser and the underside of the specimen slide. The immersion darkfield condenser has internal mirrored surfaces and passes rays of great obliquity and free of chromatic aberration, producing the best results and blackest background. Perhaps the most widely used darkfield condenser is the paraboloid, consisting of a solid piece of glass ground very accurately into the shape of a paraboloid.

As discussed above, the dry darkfield condenser is useful for objectives with numerical apertures below 0.75, while the paraboloid and cardioid immersion condensers (Figure 3) can be used with objectives of very high numerical aperture (up to 1.4). Objectives with a numerical aperture above 1.2 will require some reduction of their working aperture since their maximum numerical aperture may exceed the numerical aperture of the condenser, thus allowing direct light to enter the objective.

For this reason, many high numerical aperture objectives designed for use with darkfield as well as brightfield illumination are made with a built-in adjustable iris diaphragm that acts as an aperture stop.

This reduction in numerical aperture also limits the resolving power of the objective as well as the intensity of light in the image. Specialized objectives designed exclusively for darkfield work are produced with a maximum numerical aperture close to the lower limit of the numerical aperture of the darkfield condenser. They do not have internal iris diaphragms, however the lens mount diameters are adjusted so at least one internal lens has the optimum diameter to perform as an aperture stop.

The cardioid condenser is very sensitive to alignment and must be carefully positioned to take advantage of the very sharp cone of illumination, making it the most difficult darkfield condenser to use. In addition, the condenser produces a significant amount of glare, even from the most minute dust particles, and the short focal length may result in poor illumination on objects that exceed a few microns in size or thickness. When choosing microscope slides for quantitative high-magnification darkfield microscopy, make certain to select slides made from a glass mixture that is free of fluorescent impurities.

Careful attention should be paid to the details of oiling a high numerical aperture condenser to the bottom of the specimen slide. It is very difficult to avoid introduction of tiny air bubbles into the area between the condenser top lens and the bottom of the microscope slide, and this technique should be practiced to perfection. Air bubbles will cause image flare and distortion, leading to a loss of contrast and overall image degradation.

Problems are also encountered when using microscope slides that are either too thick or too thin. Many darkfield condensers contain the range of usable slide thickness inscribed directly on the condenser mount. If the slide is too thick, it is often difficult to focus the condenser without resorting to a higher viscosity immersion oil. On the other hand, slides that are too thin have a tendency to break the oil bond between the condenser and the slide. It is a good idea to purchase precision microscope slides of the correct thickness to avoid any of the problems mentioned above.

High numerical aperture condensers, whether intended for use dry or with oil, must be accurately centered in the optical path of the microscope to realize optimum performance.

To achieve this, many darkfield condensers are built with a small circle engraved onto the upper surface to aid in centering the condenser. Centering is performed with a low power (10x-20x) objective by imaging the engraved circle and using the condenser centering screws to ensure the circle (and condenser) are correctly centered in the optical path.


4.0 B-500DK DARKFIELD MICROSCOPE CONFIGURATION

- In a manner similar to that with low-power darkfield microscope configuration, place a specimen on the stage and focus it using the 10x objective (refer to B-500 general manual for the description of the standard bright field technique).
- Use the condenser rack knob to lower and remove the brightfield substage condenser.

- Insert the reflected light high numerical aperture condenser into the holder, and secure it with the lock screw. Slowly raise the condenser until it is less than 2 millimeters from the underside of the specimen slide.

- Select a darkfield specimen and place it onto the microscope stage between the objective and the condenser. The condenser will project a spot of light onto the sample that can be used for centering the optical path. Use the condenser centering screws to move the ring of light into the center of the viewfield. It could be useful to vary the height of the condenser in order to view the spot.


Figure 5 Centering darkfield condenser


- It is often advantageous to use a low power 10x objective when centering high numerical aperture darkfield condensers. When viewing a specimen with the 10x objective while slowly raising and lowering the condenser, a point will be reached where a bright spot will appear in the field of view as illustrated in Figure 5(a). As the condenser is slightly raised or lowered, a dark spot similar to the one shown in Figure 5(b) will appear, if the condenser is properly centered. In cases where the condenser is not properly aligned and centered, a typical field of view might look like that shown in Figure 5(c). The ideal and correct positioning of the condenser is illustrated in Figure 5(a), and the condenser should be adjusted until the field of view appears in this manner, with the condenser centering screws.
- To apply oil to the immersion condenser, remove the microscope slide and slowly rack the condenser down below the mechanical stage and place a drop of oil directly on the front lens.

- Replace the slide, then slowly raise the condenser until contact is made between the oil droplet and the bottom surface of the slide. Carefully observe the oil contact area to determine if any air bubbles have been introduced into the space between the condenser top lens and the microscope slide. It should be noted that no matter what the magnification of the objective is, if an oil immersion condenser is used without oil, light will not emerge from the condenser. Air bubbles will scatter light and should be observable by removing the eyepiece and viewing the back focal plane of the objective with a centering telescope. If no bubbles are present, proceed on to the next step. When bubbles are detected, remove all traces of oil from both the slide and the condenser top lens and repeat the procedure. Make certain both surfaces are clean and completely free of oil, dust, dirt and other contaminating artifacts before re-oiling the condenser.

- Swing the provided special 100X immersion objective to a position in the nosepiece where it will be the next objective rotated into the optical pathway directly above the specimen.
- Position the area to be imaged in the center of the microscope's optical path using lower power objective (e.g. 10x or 40x). Next, the oil immersion objective (without oil) is placed into the optical path and the specimen area of interest is roughly brought into center of the viewfield and focused. This pre-positions all of the components of the system in preparation for the addition of oil. Swing the immersion objective to an adjacent stop on the microscope nosepiece and apply oil both to the objective front lens and specimen as described in the next point.
- Use a piece of lint-free paper or an eyedropper to place a small droplet of oil directly on the front lens of the immersion objective. Next, place another droplet of oil on the cover glass immediately above the area to be imaged.
- Next, quickly rotate the oiled objective into position above the specimen merging the two oil droplets (one on the sample and one on the objective) into a single pool. A bright flash should occur (from scattered light) the instant oil on the front lens of the objective comes into contact with the microscope slide.
- Check for air bubbles by observing the back focal plane of the objective with the provided centering telescope. If bubbles are present, swing the objective to the next detent stop on the microscope nosepiece and wipe off excess oil with a piece of lint-free lens tissue. Re-apply oil to the objective front lens and try again.
- With the centering telescope, and the 100x objective inserted, check that the condenser is well centered (a symmetrical ring of light must be seen). Finely center using the condenser centering screws.

OIL
(SLIDE-OBJECTIVE)


OIL
(CONDENSER - SLIDE)

- The 100X objective has an internal iris diaphragm designed to allow adjustment of the numerical aperture. It should be adjusted at this time.


ROTATE THIS TO
ADJUST THE
INTERNAL IRIS

Remove the eyepiece and place the centering telescope into the eye tube to observe the rear focal plane of the objective. If a bright ring of light is seen around a central dark region, the aperture of the objective is too great and must be reduced using the internal iris diaphragm.

VIEW INSIDE THE TELESCOPE
WITH A WELL CENTERED
CONDENSER

Adjust the objective's iris diaphragm until the ring of light falls out of the viewfield.

VIEW INSIDE THE TELESCOPE
WHEN THE LIGHT RING IS DARKE-
NED BY CLOSING THE OBJECTI-
VE'S IRIS

If a crescent-shaped arc of light is observed, then the condenser **is not properly aligned** in center of the microscope optical axis. Use the LED centering screws to fine align the condenser so that the arc disappears or becomes a complete circle of light.

- Replace the eyepiece and focus the specimen, which should now appear brightly and evenly illuminated upon a dark background.

Careful attention should always be given to microscope alignment, irrespective of whether the illumination mode is brightfield, darkfield, phase contrast or some other contrast enhancement technique. Time spent in this endeavor will be repaid in excellent performance of the microscope both for routine observation and critical photomicrography.

There are numerous common problems associated with darkfield microscopy and photomicrography. These range from insufficient illumination and condenser misalignment to using a field stop of incorrect size. Most darkfield illumination problems are associated with the substage condenser, and this should be the first suspect when things do not work properly.

The following problems and solutions should be used as a guide when imaging specimens using this technique.

Problem There is insufficient illumination to make the specimen visible, or the specimen is visible but very faint.

Solution Check the substage condenser to ensure that it is positioned correctly. Also, check the light path for color and/or neutral density filters, polarizers, retardation plates, or any other components that might reduce illumination intensity.

Check to make certain that the condenser top lens is correctly oiled to the bottom of the microscope slide. If the immersion contact has been broken, re-apply the oil and rack up the condenser until the top lens is completely immersed in oil and in contact with the bottom of the slide.

Problem The viewfield has a dark spot in the center reducing illumination intensity, but objects in the periphery are well-illuminated and appear normal.

Solution The substage condenser is incorrectly positioned. Carefully rack the condenser up and down while observing the specimen in the viewfield. Check to make certain that both the condenser aperture and field diaphragms are set to the wide open position.

Occasionally this problem arises when microscope slides of excessive thickness (greater than the standard one millimeter) are used. Use a micrometer or a set of dial calipers to determine the thickness of the slide and the slide/cover slip combination.

Problem The image has a brightfield periphery with or without dark regions in the center.

Solution This problem is typical when the substage condenser is incorrectly centered. Revert to brightfield and re-establish the conditions of Köhler illumination, then repeat the darkfield procedure, making sure the condenser is centered.

If there is an internal iris diaphragm in the objective, reduce the iris opening size to determine if too much oblique illumination is entering the front lens of the objective.

Problem The periphery of the viewfield is bright on one side only.

Solution Check the nosepiece to ensure the objective is correctly positioned in the optical axis of the microscope. There is a detent stop that should produce a positive “click” when the objective is properly positioned.

Problem Images of the specimen are not clear and lacking in sufficient contrast and detail.

Solution The specimen might not be suitable for darkfield microscopy. Many stained specimens do not scatter enough light to be clearly imaged under darkfield conditions.

Specimen thickness may also be a problem, because very thin specimens are often not imaged well with the oblique light rays emitted from darkfield condensers. Change to brightfield or phase contrast to determine if this improves specimen contrast.

Problem Bright areas that are out of focus obstruct viewing and/or photomicrography of darkfield specimens.

Solution This is probably due to dust, hair, fibers, and/or dirt contamination of an optical surface somewhere above the condenser. Thoroughly clean the specimen slide with optical-grade tissue or cotton. Occasionally, this problem arises due to contamination on the objective front lens. Carefully loosen the objective from its seat in the nosepiece and slowly turn while observing the specimen through the eyepieces. If the obstruction rotates along with the objective, then it is probably due to dust on the front lens. Remove the objective and moisten the front lens by gently exhaling on it, then clean with lens tissue or a cotton swab wrapped on a long wooden rod.

Replace the objective and check to make certain the obstruction has been completely removed.

Problem The image of the lamp condenser partially fills the viewfield, obscuring details of the specimen.

Solution The slide may be too thin or the condenser may be adjusted too high. Measure the thickness of the microscope slide and check to make certain the oil contact is not broken. If difficulty is encountered keeping the bottom of the slide oiled, try using a thicker slide.

Problem When viewing aquatic organisms in a aqueous mount, specimens drift constantly in a single direction, making observation and photomicrography difficult.

Solution Convection currents are being created due to slow evaporation at the edges of the coverslip. Place a bead of petroleum jelly or Cytoseal mounting medium around the periphery of the coverslip to seal it firmly onto the microscope slide. Commercial products are available that can be added to the water to slow the movement of these minute organisms, making them easier to photograph.

Problem When viewing the specimen through a 10x objective, there is a oblong spot of light in the center of the viewfield.

Solution The lamp condenser may be incorrectly focused. Re-focus the condenser.

Problem The viewfield has bright spots and there are problems with achieving sharp focus of the specimen.

Solution Small air bubbles may be trapped in the oil between the top of the condenser and the bottom of the microscope slide. Remove the eyepiece to observe the back focal plane of the objective, where the bubbles will be visible. To remedy, break the oil contact and remove any residual bubbles with a cotton swab. Clean excess oil from the condenser front lens and the microscope slide, then apply a clean drop of oil to the lens. If the problem persists, thoroughly clean all optical surfaces before re-applying the oil.

Problem At high magnifications, colloidal particles display incomplete Newton rings and appear to be unevenly illuminated.

Solution The condenser is probably off-center. Change to a lower power objective and re-center the condenser.

The first rule in microscopy is to keep the optical elements completely free of dust, dirt, oil, solvents, and any other contaminants. The microscope should be kept in a low vibration smoke-free room that is clean as possible and has minimal disturbance of the circulated air. Use a dust cover on the microscope when not in use and keep all accessories in air-tight containers. Avoid using corrosive solvents to clean any part of the microscope, and use only diluted soapy water to clean non-optical surfaces. Oil should be used keeping the following techniques and precautions in mind:

- The first step is to locate the sample with a low-power objective and position the microscope slide so that the area of interest is squarely in the central portion of the stage opening. It is often desirable to rotate the immersion objective (without oil) into place to ensure the correct location of the specimen. To achieve optimum results when using the oil immersion technique it is important to also apply oil to the underside of the microscope slide and to the top lens of the substage condenser to form an oil bead between the two, as described below.

For applying oil to the front lens of the substage condenser, use the following considerations:

- Make certain that the mechanical stage of the microscope has an opening of sufficient size to allow scanning of the oiled specimen without spilling oil onto the bottom of the stage. Many advanced research microscopes have a stage insert that is removable for use with oil. Remove the specimen slide from the microscope stage before applying oil to the condenser top lens.
- In a manner similar to oiling an objective, place a drop of oil on the bottom of the microscope slide and apply a similar drop to the front lens of the substage condenser. Make sure that the condenser is slightly lowered in its rack before replacing the microscope slide on the stage. After the slide has been placed on the stage, rack the condenser back to the proper position while carefully observing the oil droplets on the condenser lens and microscope slide. These droplets should merge into one as the condenser is brought into proper position.
- Swing the immersion objective to a position in the nosepiece where it will be the next objective rotated into the optical pathway directly above the specimen.
- Position the area to be imaged in the center of the microscope's optical path using lower power objectives. Next, the oil immersion objective (without oil) is placed into the optical path and the specimen area of interest is roughly brought into center of the viewfield and focused. This pre-positions all of the components of the system in preparation for the addition of oil. Swing the immersion objective to an adjacent stop on the microscope nosepiece and apply oil both to the objective front lens and specimen as described in the next point.

- Use a piece of lint-free paper or an eyedropper to place a small droplet of oil directly on the front lens of the immersion objective. Next, place another droplet of oil on the cover glass immediately above the area to be imaged.
- Next, quickly rotate the oiled objective into position above the specimen merging the two oil droplets (one on the sample and one on the objective) into a single pool. The specimen is now ready for observation and photomicrography.
- As you scan the surface of an oiled microscope slide, periodically observe the position of the stage with respect to the objective and condenser front lens to assure that there is a continuous bead of oil between these elements and the microscope slide. When scanning over a large area on a microscope slide, it is often necessary to reapply some oil to assure continuity of the imaging medium.

Keep The Lenses Clean and Oil-Free

Always strive to make sure that your microscope objective and substage condenser front lenses are kept clean and free of immersion oil. Use a Q-tip or lint-free paper or cloth to gently wipe excess oil from the lenses and microscope slide after experiments are finished. After the oil is removed, use a suitable solvent to remove traces of immersion oil from the lenses. Failure to clean lenses properly may result in small crystallites forming on the coated surfaces when immersion oils either dry or collect dust and other contaminating particles from the atmosphere.

6.0 ELECTRICS

Mains Power:	90-240 V, 50/60 Hz
Light source:	3W high-efficiency white LED

Selecting the power supply

- On the rear of the microscope you find a 3-positions switch to control the power supply of the illuminator.

Position I: power supply from line voltage


Position 0: illuminator OFF

Position II: power supply from internal battery

- A full charge of the battery provides about **2 hours** of continuous use at maximum light intensity, or about **3 hours** at medium intensity.
- The typical lifetime of the internal battery pack is about **500 charge/discharge cycles**.

Recharging the battery

- In order to guarantee an optimum battery life and efficiency, always follow these steps to recharge:
 - ◆ Connect the charger to 230V line voltage (WARNING: not compatible with 110Vac)
 - ◆ Check the settings for voltage and charging current:

Voltage: **8,4V**

Current: **1000mA**

NEVER change these settings in order to prevent damages to the battery pack.

- ◆ Connect the charger to the jack at the rear:

- ◆ Press the “discharge” button until the charger’s LED becomes yellow.
- ◆ After the discharging phase, the charger will automatically start the charging process and the LED will become red.
- ◆ When the batteries will be fully charged, the LED will become green (this could take about **3 hours**). The charger will apply a maintenance charge as long as it will be connected to the microscope.
- ◆ Disconnect the charger’s jack from the rear.
- ◆ Disconnect the charger from 230V line voltage.

Please note:

- After 2 hours of use at full intensity, stop using the battery power supply even if the LED illuminator is still on: excessive discharge may damage the batteries.
- Batteries are subject to self-discharge: after 30 days of inactivity a full charge drops to 40% of capacity, so it’s suggested to recharge the battery pack as described above.
- While the battery pack is charging, it’s still possible to use the microscope from line voltage, setting the main switch to position I.

Replacing the internal battery pack



Use ONLY this type of battery: AA NiMH 1.2V, capacity at least 2500mAh

After about 500 cycles of charge/discharge, the internal batteries lose their capacity and the autonomy begins to decrease. In order to replace the batteries, follow these steps:


7.0 USING THE RECHARGEABLE BATTERIES

- Disconnect the power plug from the microscope rear.
- Unscrew the four screws that fix the back plastic cover:


- Remove the plastic cover, and gently lay it on the back of the microscope (don't pull it too hard since it has cables attached to it).
- Inside at the bottom you can see the battery pack holder:

- Pull the battery holder toward yourself, in order to extract it.
- Remove the two-poles connector from the battery-holder:

- Remove the 6 batteries inside the holder and replace them with the new one (**pay attention to the correct polarity of each battery**).
- Reconnect the two-poles connector and replace the battery holder inside the microscope.
- Close the back plastic cover with the four screws.


Please note:

Use ONLY this type of battery: AA NiMH 1.2V, capacity at least 2500mAh.

Respect your local regulation for the disposal of nickel-metal hydride batteries.

The appliance reports the symbol:

This symbol means that the appliance can be a precious source of raw materials. Therefore, it must not be disposed of as waste, but separately collected for the recycling and the recovery of the materials it contained in it. Such materials, if improperly dispersed into the environment, can be harmful to the environment and to human health.

The producer of the equipment, Optika Microscopes, recovers, re-uses and recycles the raw materials contained in the equipment. Such recovery, however, needs your help.

When, at the end of its operating life, you shall decide to dispose of the apparatus, do not try to open it, nor to use parts of it in ways other than reported in this Manual, but bring it back to the Distributor you bought it from, or to other Optika Microscopes distributors. The Distributor shall collect the apparatus free of charge.

The recovery of the raw materials shall then take place in accordance with the European Directive 2002 / 96 / EC and all other relevant Directives. Never disassemble, nor dispose of as waste, apparatus reporting the “crossed bin” symbol indicated above.

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