Modular Series

INSTRUCTION MANUAL

Model
CLD-01

v 1.0 2018

SUMMARY

Warning

Intended use

Overview

Assembly

Using the illuminator

Equipment disposal

Warning

This cold light illuminator is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use.


We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users.

We decline any responsibility deriving from incorrect instrument use uses that does not comply with this manual.

Intended use

For research and teaching use only. Not intended for any animal or human therapeutic or diagnostic use.

Overview (with CL-11.1 sold separately)

5. Assembly of optional accessories (CL-11.1 or CL-12)

1. Insert the light guide connector in the light exit of the illuminator and lock the screws (Fig. 1-3)
2. Insert the focus lens on top of the light guide arms (Fig. 4)
3. Connect the power supply. (Fig. 5)

4. Connecting polarizing filters (Fig. 6)
Insert the polarizing filter on top of the focussing lens and lock the fixing screw.
Adjust the polarizing effect rotating each filter according to the analyzer installed on the microscope head.

6. Using the Illuminator

Turn on the illuminator moving the light intensity knob ① (Fig. 7).

Set the desired light intensity according the specimen under examination.


Bend the self-supporting arms ② to aim the light beam (Fig. 8).

Each arm can be oriented in a separate way, in order to achieve a proper illumination (Fig. 9)

Equipment disposal

Art.13 Dlsg 25 july 2005 N°151. "According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal."

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste.

The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection.

The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment.

Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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OPTIKA® S.r.l.