



LCD Mikroskop 50x -2000x / LCD Microscope 50x -2000x

- Revolverkopf mit 3 hochwertigen Objektiven(4x / 10x / 40x)
- LED-Beleuchtung für Auflichtmikroskopie (230V Stromanschluss)
- LED-Beleuchtung für Durchlichtmikroskopie (230V Stromanschluss)
- Beleuchtungsdimmer und Schalter für Auf- und Durchlicht
- 6-Farben-Filterscheibe zur optimalen Betrachtung von farblosen Präparaten
- Präzisions-Kreuztisch für genaue Ausrichtung des Beobachtungsobjekts
- Präparate (Garneleneier, Seewassersalz und Hefe) / Fixiermittel
- Mikroskopiebesteck / Garnelenbrutanlage / Mikrotom (zum Anfertigen von Schnitten)
- Box mit Präparaten, leeren Objektträgern und Deckgläsern

- Revolver head with three high-quality objectives (4x / 10x / 40x)
- LED illumination for Stereo microscopy (230V)
- LED illumination for transmission-light-microscopy (230V)
- Switch with dimmer for both light sources
- 6 colour filter slides for optimal viewing of colourless preparations
- · Precision mechanical desk for accurate orientation of the object
- · Preparations (shrimp eggs, yeast, sea salt) / Gum media
- Microscope accessory kit / hatchery / microtome (for thin preparation cuts)
- Box with preparations, blank objective slides and cover glasses