



* Head not included

Model:
SZ-STL1H

Typology:
STAND FOR MODULAR STEREO MICROSCOPES

Description:
Simple overhanging stand with head holder and focusing mechanism (head sold separately).

Illumination	This stand has no illumination system, but it can be used with external illuminators (halogen, LED, optical fibers).
Focusing	By means of a rack and pinion controlled by a pair of knobs placed on both sides of the head holder. Tension of focusing can be adjusted. The head holder can slide along a pillar. The stereo head connection is 76mm.
Dimensions	Simple overhanging stand with pillar height 42 cm, diameter 32 mm. Horizontal arm: 43 cm. Base: 23x23x8h cm. Weight: 15 Kg.