


## LCACHN 20XIPC UPLFLN-IPC / CACHN-IPC / LCACHN-IPC




This series is used in combination with the pre-centered phase contrast slider CKX3-SLP. When changing the objective magnification in phase contrast observation, no centering adjustment is necessary.

## **Specifications**

opecinications	
Numerical Aperture (NA)	0.4
Working distance (W.D.) (mm)	3.2
Field Number (FN)	22
Cover Glass Thickness (mm)	1
Immersion	-
Spring	-
Correction Collar	-
Iris Diaphragm	-
Waterproof & Oil proof function	-
Ultra Wide anti-reflection coating	-
Super apochromatic	-
Brightfield	Good
Darkfield	Good
DIC	N/A
Phase Contrast	Excellent
Relief Contrast	N/A
Polarized Light	Limitation
Fluorescence (B, G Excitation)	Good
UV Fluorescence (at 365nm)	N/A
IR DIC	N/A
Multi photon excitation	N/A
	·

## **LCACHN 20XIPC**

## **Dimension**



• Specifications and appearances are subject to change without any notice or obligation on the part of the manufacturer.

