



BA210 Binocular LED

The BA210 is designed for both educational and teaching environments from basic life sciences to medical applications. Through the Motic Infinity Optics (CCIS®), the BA210 delivers a new higher level of performance in education and training. Contrast methods like Phase Contrast, Polarization and Darkfield are easily performed by using optional accessories.

Model BA210 Binocular LED

Optical System Colour Corrected Infinity Optical System (CCIS®)

Observation tube Binocular head, Siedentopf type

Interpupillary distance 55-75mm

Inclination 30° inclined, 360° rotating.

Eyepieces N-WF10X/20mm with diopter adjustment, +/- 5 diopter

Nosepiece Reversed quadruple revolving nosepiece

Objective classification Infinity Corrected CCIS EC Plan Achromatic DIN 45mm

Objectives 4X/0.10 (WD 15.9mm), 10X/0.25 (WD 17.4mm), 40X/0.65/S (WD 0.5mm),

100X/1.25/S-Oil (WD 0.15mm)

Objective mounting thread W 4/5" x 1/36" (RMS standard)

Stage Built-in low position coaxial mechnical stage with sample holder

Stage size 140X135mm

Mechanical stage X&Y range 76X50mm

Upper limit stop Upper limit stop preset but adjustable

Condenser Focusable Abbe Condenser N.A. 1.25 with iris diaphragm and slot for contrast sliders

Focus mechanism Coaxial coarse and fine focusing system with tension adjustment

Fine Focus precision 2 µm minimum increment

Z-axis movement 25mm

Filter holder On top of the illuminator

Illumination 3W LED illumination with intensity control

Transformer Internal

Power supply 100-240V (CE)

Accessories included Immersion oil 5ml, power cord, cord hanger, dustcover, allen key, thumb screw and spare fuse

Dimensions 360x190x395mm

Weight 7,2kg



Add a Moticam camera to capture, document, annotate and share images and videos with the Motic Images Plus software.