## KERN

## Compound microscope KERN OBF-1





Trinocular version



Simple polarising attachment

## LAB LINE

# The high-performance compound microscope for every laboratory with fixed, pre-centred Koehler illumination

### **Features**

- The KERN OBF models are excellent, stable laboratory microscopes for all common routine applications. A central feature of this adaptable, robust microscope series is the stable mechanism which can be adjusted precisely
- Depending on the application, there is a choice of models with strong, continuously dimmable 3 W LED or 20 W halogen illumination (Philips)
- The fixed, pre-centred and focusable 1,25
   Abbe condenser with aperture diaphragm and field diaphragm gives you a simplified Koehler illumination, without having to move the centre
- The large mechanical stage and its specimen holder holds up to two samples at the same time and is quick and easy to focus using a coaxial coarse and fine focusing knob on both sides

- A large selection of eyepieces, objectives and colour filters as well as a darkfield condenser and a simple polarising unit are available to you as accessories
- A protective dust cover, eye cups, as well as multi-lingual user instructions are included in the scope of delivery
- A C-mount adapter is required to connect a camera to the trinocular version. You can select this adapter from the following model outfit list
- Please find detailed information in the following model outfit list

## Scope of application

 Haematology, urology, gynaecology, dermatology, pathology, microbiology and parasitology, immunology, oncology, entomology, vets, water analysis and breweries

## Applications/Samples

 Translucent, thin, low-contrast, challenging samples (e.g. living mammal cells, bacteria, tissue)

## Technical data

- · Finite optical system (DIN)
- Quadplex nosepiece
- Siedentopf 30° inclined/360° rotatable
- Diopter adjustment: One-sided
- Overall dimensions W×D×H 395×200×380 mm
- Net weight approx. 6,7 kg

STANDAR	D							OPTION		
Ø	00		Δ	Ф	Ð	<b>-</b>		•		<
360°	BINO	TRIN0	ABBE	HAL	LED	230 V	1 DAY	DF	POLAR	S

Model		Standard configuration				
KERN	Tube Eyepiece		Objective quality	Objectives	Illumination	
OBF 121	Binocular	HWF 10×/Ø 18 mm	Achromatic		20 W Halogen (transmitted)	
OBF 122	Binocular	HWF 10×/Ø 18 mm	Plan		20 W Halogen (transmitted)	
OBF 123	Binocular	HWF 10×/Ø 18 mm	Plan	4×/10×/40×/100×	3 W LED (transmitted)	
OBF 131	Trinocular	HWF 10×/Ø 18 mm	Achromatic	4^/10^/40^/100^	20 W Halogen (transmitted)	
OBF 132	Trinocular	HWF 10×/Ø 18 mm	Plan		20 W Halogen (transmitted)	
OBF 133	Trinocular	HWF 10×/Ø 18 mm	Plan		3 W LED (transmitted)	

# **KERN OPTICS CATALOGUE 2020**



## Compound microscope KERN OBF-1

Model outfit				Model	KERN		Order number		
			OBF OBF OE				-		
	HWF 10×/ø 18 mm	121	131	122 <b>√√</b>	132	123	133	OBB-A1403	
	WF 16×/Ø 13 mm	00	00	00	00	00	00	OBB-A1463	
Eyepieces (23,2 mm)	HWF 10×/Ø 18 mm (with Pointer)	0	0	0	0	0	0	OBB-A1334 OBB-A1348	
	HWF 10×/Ø 18 mm (reticule 0,1 mm) (non-adjustable)	0	0	0	0	0	0	OBB-A1349	
	4×/0,10 W.D. 18,6 mm	<b>✓</b>	<b>✓</b>					OBB-A1111	
		-							
	10×/0,25 W.D. 6,5 mm	<b>*</b>	<b>√</b>					OBB-A1118	
Achromatic objectives	40×/0,65 (spring-loaded) W.D. 0,47 mm	<b>*</b>	<b>√</b>					OBB-A1112	
,	100×/1,25 (oil) (spring-loaded) W.D. 0,07 mm	<b>✓</b>	<b>✓</b>					OBB-A1109	
	20×/0,40 (spring-loaded) W.D. 1,75 mm	0	0					OBB-A1110	
	60×/0,85 (spring-loaded) W.D. 0,1 mm	0	0					OBB-A1113	
	4×/0,10 W.D. 14,5 mm			✓	✓	✓	✓	OBB-A1255	
	10×/0,25 W.D. 5,65 mm			✓	✓	✓	✓	OBB-A1238	
	40×/0,65 (spring-loaded) W.D. 0,85 mm			<b>✓</b>	✓	<b>✓</b>	✓	OBB-A1256	
Plan objectives	100×/1,25 (oil) (spring-loaded) W.D. 0,07 mm			<b>✓</b>	✓	✓	✓	OBB-A1239	
	20×/0,40 (spring-loaded) W.D. 1,5 mm			0	0	0	0	OBB-A1249	
	60×/0,85 (spring-loaded) W.D. 0,07 mm			0	0	0	0	OBB-A1269	
	100×/1,0 (water) (spring-loaded) W.D. 0,18 mm	0	0	0	0	0	0	OBB-A1441	
Binocular tube	<ul> <li>Siedentopf 30° inclined/360° rotatable</li> <li>Interpupillary distance 50 - 75 mm (for non-infinity system)</li> <li>Diopter adjustment: One-sided</li> </ul>	~	0	<b>✓</b>	0	<b>✓</b>	0	OBB-A1129	
Trinocular tube	Siedentopf 30° inclined/360° rotatable     Interpupillary distance 50 – 75 mm     Light distribution 20:80 (for non-infinity system)     Diopter adjustment: One-sided	0	~	0	~	0	~	OBB-A1345	
Mechanical stage	Stage size W×D 145×130 mm     Travel 76×52 mm     Coaxial coarse and fine focusing knobs, scale: 2 µm     Two slide holder	~	<b>✓</b>	<b>✓</b>	<b>✓</b>	<b>*</b>	<b>✓</b>		
Condenser	Abbe N.A. 1,25 precentered (aperture diaphragm)	1	✓	1	✓	✓	✓	OBB-A1103	
Darkfield condenser	N.A. 0,85 – 0,91 (dry, paraboloid)	0	0	0	0	0	0	OBB-A1422	
	20 W Halogen spare bulb (transmitted)	1	1	1	1			OBB-A1370	
Illumination	3 W LED illumination system (transmitted) (non-rechargeable)					~	<b>✓</b>		
Polarising unit	Analyser/Polariser	0	0	0	0	0	0	OBB-A1277	
	Blue (built-in)	1	✓	✓	✓	✓	✓		
Colour filters	Green	0	0	0	0	0	0	OBB-A1188	
for transmitted illumination	Yellow	0	0	0	0	0	0	OBB-A1165	
	Grey	0	0	0	0	0	0	OBB-A1183	
	0,47× (focus adjustable)				0		0	OBB-A1135	
	0,5× (focus adjustable)		0					OBB-A1515	
C-Mount					0		0	OBB-A1142	
	1×		0					OBB-A1514	

✓ = Included with delivery

O = Option

## **KERN OPTICS CATALOGUE 2020**

## **Pictograms**



360° rotatable microscope head



Fluorescence illumination for compound microscopes With 3 W LED illumination and filter



#### WLAN data interface

For transmitting of the picture to a mobile display device



Monocular Microscope

For the inspection with one eye

For the inspection with both eyes



Phase contrast unit For a higher contrast



**HDMI** digital camera

For direct transmitting of the picture to a display device



Binocular Microscope



Darkfield condenser/unit

For a higher contrast due to indirect illumination



PC software

To transfer the measurements from the device to a PC.



Trinocular Microscope

For the inspection with both eyes and the additional option for the connection of a camera



Polarising unit

To polarise the light



Automatic temperature compesation

For measurements between 10 °C and 30 °C



Ф

Abbe Condenser

Halogen illumination

With high numerical aperture for the concentration and the focusing of light

For pictures bright and rich in contrast



Infinity system

Infinity corrected optical system



Protection against dust and water

splashes IPxx



The type of protection is shown by the pictogram.





Zoom magnification

For stereomicroscopes

Parallel optical system For stereomicroscopes, enables



**Battery operation** 

Ready for battery operation. The battery type is specified for each device.



long-life illumination Incident illumination

LED illumination

For non-transparent objects



fatigue-proof working

Integrated scale In the eyepiece



**Battery operation rechargable** Prepared for a rechargable battery

operation





SCALE

SD card For data storage



Mains adapter

230V/50Hz in standard version for EU. On request GB, AUS or USA version.



Fluorescence illumination

For stereomicroscopes



USB 2.0 digital camera

For direct transmitting of the picture to a PC



Power supply

Integrated in microscope. 230V/50Hz standard EU. More standards e.g. GB, AUS or USA on request.



Fluorescence illumination for compound microscopes

With 100W mercury lamp and filter



USB 3.0 digital camera

For direct transmitting of the picture to a PC



Package shipment

The time required to manufacture the product internally is shown in days in the pictogram.

## **Abbreviations**

Adapter for the connection of a C-Mount

camera to a trinocular microscope

LWD Long Working Distance SWF Super Wide Field (Field number at

least Ø 23 mm for 10× eyepiece)

**FPS** Frames per second N.A. Numerical Aperture

Working Distance W.D.

H(S)WF High (Super) Wide Field (Eyepiece with high eye point for wearers of glasses)

SLR

Single-Lens Reflex camera

WF Wide Field (Field number up to

Ø 22 mm for 10× eyepiece)

## Your KERN specialist dealer: