



*Simple stand, pillar type*

<b>Observation Method - Incident Light</b>	Brightfield	As optional
	Simple polarized light	As optional
<b>Stand</b>	Type	Pillar
	Pillar diameter (mm)	32
<b>Stage</b>	Black/white contrast plate for incident ill.	Yes, diameter 95 mm
	Specimen clips	Yes
<b>Accessories Included</b>	Dust cover	Yes
	User Manual	Digital version (downloadable)
<b>Product Dimensions</b>	Total height (mm)	275
	Base height (mm)	25
	Base width (mm)	210
	Base depth (mm)	270
<b>Product Weight</b>	(kg)	1.4