



Base with oblique illumination

Observation Method - Transmitted Light	Brightfield	Yes
	Diascopic illumination	Yes, with rotating mirror and diffuser
Stand	Type	Pillar
	Pillar diameter (mm)	32
Stage	Frosted diffusion plate for transmitted ill.	Yes, diameter 100 mm
	Specimen clips	Yes
Transmitted Illumination	Type	X-LED
	X-LED type	X-LED3
	Light source power (W)	3.6
	Illuminance (lux)	8,000
	Brightness control	Manual
	Lifetime (hours)	> 65,000
	Temperature (K)	6,300
	Max. required power (W)	5
Power Supply for Illumination	Type	External
	Microscope connector	Jack, 2.1 mm
	Power plug type	Multi-plug (EU, UK, US)
	Input voltage	100/240 Vac, 50/60 Hz
	Output voltage	6 Vdc 2.5 A
Accessories Included	Dust cover	Yes
	User Manual	Digital version (downloadable)
Product Dimensions	Total height (mm)	340
	Base height (mm)	85
	Base width (mm)	210
	Base depth (mm)	270
Product Weight	(kg)	2,4